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## Bibliografia

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# rubagen

## Rapid latex particle agglutination test on slide for the qualitative and quantitative determination of rubella virus antibodies in serum.

### Summary

Rubella is a mild viral disease affecting children and adults. However, when contracted in the first trimester of pregnancy, rubella virus may infect the fetus through the placenta causing multiple abnormalities commonly referred to as "congenital rubella syndrome". Other consequences of rubella infection during pregnancy may include spontaneous abortion and stillbirth (1-3).

The availability of an attenuated rubella virus vaccine has greatly reduced the natural incidence of rubella infection. Nevertheless, it is recommended that all women of childbearing age be tested for the presence of rubella antibody to assure that nonimmune individuals are detected and subsequently vaccinated (4-5).

The detection of rubella immune status was performed during many years by various serological methods, specially the hemagglutination inhibition test (HAI) (6-8). However, latex agglutination in comparison with HAI is quicker and easier to perform. The sensitivity of **rubagen** correlates perfectly with that of the HAI test.

### Principle

The **rubagen** reagent is a suspension of polystyrene latex particles of uniform size coated with soluble rubella inactivated virus antigen. Latex particles allow visual observation of the antigen-antibody reaction. If rubella virus antibodies are present in the serum sample, the latex suspension changes its uniform appearance and a clear agglutination becomes evident.

When the latex reagent is mixed with the serum, if the serum contains approximately more than 10 IU/ml of rubella virus antibodies, a clear agglutination will appear. Results are expressed in IU/ml of total anti-rubella virus antibodies based on the WHO International Standard (1970).

### Components

**REF** 3000-4001.

- a) **REAG** Latex reagent: 1 x 1.7 ml.  
Suspension of polystyrene latex particles coated with rubella virus antigen in a buffer. Contains < 0.1% sodium azide.
- b) **CONTROL +** Positive control: 1 x 1.0 ml.  
Diluted positive human serum. Ready to use. Contains < 0.1% sodium azide.
- c) **CONTROL -** Negative control: 1 x 1.0 ml.  
Diluted normal human serum. Ready to use. Contains < 0.1% sodium azide.
- d) **DIL** Dilution buffer: 1 x 50 ml.  
Phosphate buffered saline, containing bovine serum albumin. Contains < 0.1% sodium azide.
- e) **SLIDES** Disposable slides: 1 x 13 units.

### Precautions

**rubagen** is intended for IN VITRO diagnostic use.

Sodium azide may react with lead or copper pipes and plumbing creating highly explosive metal azides. Flush drains with water thoroughly after disposing of the remains of reagents.

All human source material used in the preparation of this product was found to be negative for the presence of HIV-1/HIV-2 and HCV antibodies, as well as for the hepatitis B surface antigen, using a method approved by the Food and Drug Administration (USA).

WARNING: POTENTIALLY BIOHAZARDOUS MATERIAL.

Because no test method can offer complete assurance of the absence of infectious agents, this product should be handled with caution (9).

Dispose all used materials in a suitable biohazardous waste container.

#### Hazard class

Not classified.

#### Hazard statements

None.

#### Precautionary statements

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+351+338: IF IN EYES rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

 Biokit, S.A. - www.biokit.com  
Can Malé s/n - 08186 Lliçà d'Amunt - Barcelona - Spain

 **Biokit**  
A Werfen Company

  
0843

3800-0083  
R03 06/17

## Storage

The reagents will remain stable until the expiration date shown on the label, if stored between 2 and 8°C. Do not freeze.

## Material required not provided

Test tubes, automatic pipettes, rotator, stirrers and timer.

## Sample collection

Use fresh serum. Samples can be stored at 2-8°C for 8 days. For longer periods, samples should be frozen (-20°C). For diagnosis of acute rubella infection, paired sera (acute and convalescent) should be obtained. The first sample should be collected as early as possible after rash onset or exposure and the second one should be obtained 1-3 weeks later. Both sera should be tested simultaneously using the quantitative procedure (7). For qualitative antibody assay a single sample is sufficient.

## Procedure

Quality control: Before performing a set of determinations it is advisable to test the reagent with each of the controls included in the kit. Both controls should be used following the steps outlined in the QUALITATIVE TEST but without doing any dilution as they are prediluted. The reaction between the positive control and the reagent should show a clear agglutination, different from the uniform appearance of the negative control. If these results are not obtained, do not use the kit.

To ensure drops of correct volume, wipe dry the dropper tips with a tissue prior to use. For correct delivery, reagent and controls droppers must be held vertically.

## QUALITATIVE TEST

- Allow the reagents to reach room temperature.
- Dilute the serum with the dilution buffer and mix well. The following table shows the different serum dilutions to be performed depending on the sensitivity of the latex reagent which is printed on the outside of the box.

Sensitivity level of 10 IU/ml		
Latex sensitivity	Dilution	
2.5 IU/ml	1:4	(50 µl serum + 150 µl buffer)
2.0 IU/ml	1:5	(50 µl serum + 200 µl buffer)
1.5 IU/ml	1:6.7	(50 µl serum + 285 µl buffer)

- Place 25 µl of the serum dilution (or one drop of control) onto one section of the slide.
- Shake the reagent vial and add one drop of reagent next to the drop of diluted serum.
- Mix both drops with a stirrer covering the whole surface of the slide section.
- Rotate the slide for 8 minutes manually or on a rotary shaker set at 80-100 rpm.
- Observe for presence or absence of agglutination.

## Reading results

POSITIVE REACTIONS: Presence of agglutination.

- 3+ Large clumping with clear background.
- 2+ Moderate clumping with fluid slightly opaque in background.
- 1+ Small clumping with opaque fluid in background.

NEGATIVE REACTIONS:

Absence of agglutination, uniform suspension.

## Interpretation of the results

The presence of any clear agglutination indicates the presence of antibodies against rubella virus in the serum. This indicates previous exposure to the virus and the immunization of the individual against it.

The absence of agglutination indicates the absence of antibodies to rubella virus in the serum or their presence in a very low or no detectable level. In this case, the individual can be infected or reinfected by the virus.

In pregnant women, specially within the first trimester, it is recommended to repeat the test a month later in order to detect any seroconversion.

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## QUANTITATIVE TEST

- Place 25 µl of dilution buffer onto slide sections 1 through 8.
- Using an automatic pipette, place 25 µl of the serum dilution performed in the qualitative test onto slide section 1.
- Using the same pipette take in and release the serum and the dilution buffer on section 1 several times until they are well mixed.
- Take 25 µl of the mixture made on section 1 and transfer it to section 2.
- Repeat the aforementioned operations to obtain a thorough mixing of reagents, through section 8, thereafter discarding 25 µl.
- Test each dilution as described in the section QUALITATIVE TEST.

## Interpretation of the results

The approximate titer will correspond to the highest serum dilution that still presents a clearly visible agglutination. When expressing the titer in IU/ml the sensitivity of the latex reagent which is printed on the outside of the box must be multiplied by the reciprocal of the last dilution of serum giving a positive result. The sensitivity of the latex reagent has been determined against the 1st WHO International Standard.

- e.g. - Latex reagent sensitivity: 2 IU/ml.
- Last positive dilution: 1:20.
- Titer: 20 x 2 IU/ml = 40 IU/ml.

The diagnosis of acute rubella infection should be confirmed serologically by the presence of specific IgM antibody or by comparing antibody titers in paired sera. The first sample should be collected as early as possible after rash onset or exposure, while the second sample should be obtained 1-3 weeks later. A four-fold or greater rise in titer between the first and the second sample is considered diagnostic of acquired rubella infection (2-4).

Also a seroconversion may be seen after a vaccination procedure.

## Limitations of the procedure

- Test results obtained with **rubagen** must be evaluated by the physician in light of the clinical symptoms shown by the patient.
- Acute and convalescent sera must be tested at the same time. The absence of a four-fold rise in titer does not exclude the possibility of exposure and infection.
- In one evaluation, a prozone effect was observed with one strong positive IgM serum sample. If this possibility is suspected, the test should be repeated using a 1:20 dilution of the serum.

## Expected values

Antibodies against natural infection or vaccine rubella virus are prevalent in 80% of young adults. Adults are rarely afflicted with rubella, but women of childbearing age can acquire primary rubella infection by intensive contact with young children (3).

The prevalence data reported for rubella vary depending of the country. A revision study reported seroprevalence of rubella in women of childbearing age varying from 81.1% (Korea), 83.9 (Cameroon), 89.6% (Spain) to 94.3% (Switzerland) (10).

## Performance characteristics

The **rubagen** sensitivity was studied in comparison to a commercially available IgG and IgM Microparticle Enzyme Immunoassays (MEIA) using a serum panel of 77 positive samples for IgG antibody, 16 positive samples for IgG and IgM rubella antibodies, and 69 negative samples. The relative sensitivity was 100%, all 93 samples were positive with the latex reagent. The relative specificity was 76.8% (53/69). The 16 discrepant samples were tested with a quantitative commercial assay that detected 8 of them as positive, 7 as "doubtful" (6-9 IU/ml), and 1 sample as negative. Eliminating the doubtful samples from calculation, the specificity of the test was 98.1%. Once resolved the discrepancies, the positive predictive value for **rubagen** was 99% and the negative predictive value was 100%.

A comparative study between **rubagen** and two other commercially available latex agglutination tests was done using 19 serum samples classified as negative (< 10 IU/ml) or equivocal (10-15 IU/ml), 34 positive (> 15 IU/ml) and 10 IgM-positive by an EIA technique. The samples were tested following the instructions of each manufacturer. **rubagen** (10 IU/ml of detection) was the most sensitive kit detecting all the equivocal EIA samples (10/10) and 4 negative samples (7-10 IU/ml). A second latex reagent (10 IU/ml of detection) detected 9/10 of the equivocal samples and 2 negative samples (8 IU/ml). The third latex reagent (3 IU/ml of detection) detected 6/10 of equivocal samples and 2 negative samples (8 IU/ml). The 34 positive serum samples with EIA results ranging from 15 to 200 IU/ml were found positive by the 3 latex reagents, except 2 serum samples that were negative by the third latex reagent. This third reagent also did not detect 4/10 serum samples (> 200 IU/ml) because of a prozone effect. The 10 serum samples positive for specific IgM were also tested, **rubagen** and the second reagent detected 9 of the samples while one sample showed a prozone effect with both reagents. The third reagent showed prozone effect with 4 of the serum samples.

In an external evaluation of multiple kits for detection of rubella antibodies, 125 samples previously characterized were tested. All 95 negative samples tested were also negative with **rubagen**. The 35 positive serum samples, 20 of them classified as "weakly" positive (15-40 IU/ml) and 10 as "strong" positive (> 500 IU/ml), were found also positive with this kit. Five seroconversion panels of vaccinated persons, with 4 consecutive samples each obtained one before and 3 after the vaccination, were tested and **rubagen** detected correctly all of them.

In another study **rubagen** accurately detected the seroconversion of paired sera from 10 persons with natural infection and 4 who had been vaccinated. All serum pairs showed a four-fold or greater rise in titer.

#### Precision

10 serum samples were quantified 3 times on different days. The results indicate that the reproducibility of the reagent (within one dilution) was 100%.

## References

1. Gregg NM. Congenital cataract following German Measles in the mother. Trans Ophthalmol Soc Aust 3: 35-46, 1941.
2. Chernesky MA, Mahony JB. Rubella virus. In "Manual of Clinical Laboratory Immunology". Third Edition. American Society for Microbiology, 536-539, 1986.
3. Doerr HW. Viral diseases. In: Thomas L. Clinical laboratory diagnosis. Use and assessment of clinical laboratory results. 1<sup>st</sup> edition. TH-Books, Frankfurt/Main, Germany, 1998.
4. Rubella prevention - Recommendations of the immunization practices advisory committee (ACIP). November 23, 1990 / 39(RR 15); 1-18.
5. Skendzel LP. Rubella immunity. Defining the level of protective antibody. Am J Clin Pathol 106: 170-173, 1996.
6. Stewart GL, Parkman PD, Hopps HE, Douglas RD, Hamilton JP, Meyer HM. Rubella-Virus hemagglutination-inhibition test. New Eng J Med 276: 554-557, 1967.
7. Meegan JM, Evans BK, Horstmann DM. Comparison of the latex agglutination test with the hemagglutination inhibition test, enzyme-linked immunosorbent assay, and neutralization test for detection of antibodies to Rubella virus. J Clin Microbiol 16: 644-649, 1982.
8. Grangeot-Keros L, Pillot J. Étude critique du sérodiagnostic de la Rubéole: Évaluation comparative des méthodes classiques et nouvelles et leur signification immunitaire. Bull de L'Institut Pasteur 83: 375-388, 1985.
9. Biosafety in Microbiological and Biomedical Laboratories. CDC/NIH manual, 5th edition, 2007.
10. Zufferey J, Jacquier P, Chappuis S, Spinnler O, Hohlfeld P, Zuber PLF, Bille J. Seroprevalence of rubella among women of childbearing age in Switzerland. Eur J Clin Microbiol Infect Dis 14: 691-696, 1995.

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