biokit

08186 Llissà d'Amunt
Barcelona
Spain

Tel.:+ 34 93 860 90 00
Fax:+34 93 860 90 17
e-mail: biokit@biokit.com
www.biokit.com

Date: 02.02.09

Ref: 271/MKT/Q/50

SUBJECT **Quantex HbA1c. NGSP factor.**

PRODUCT EVALUATION

Quantex HbA1c evaluation with NGSP method. Implementation of NGSP factor in reagent insert and control values.

Quantex HbA1c% (REF. 3000-2314) kit was evaluated in an Ilab 600 Analyser (Instrumentation Laboratory) and compared with a method standardized against NGSP.

A multiplying factor for this test was obtained after analysing all data reported by the METHOD COMPARISON STUDIES. This factor can be used with our product in order to obtain the same standardization as the reference method used in this evaluation.

We believe that this information can be really interesting for the quantex clients who would be interested in applying the NGSP standardization. That's why we are implementing this factor in our quantex HbA1c% reagent insert.

At the same time, with the aim to facilitate the work, we have modified the control values insert sheet so that it shows the expected ranges both considering and not considering this factor.

Customers who prefer to continue with our standardization will not have to do any adjustments.

As always, we will include a note in our kits (in yellow paper) explaining these implementations.

This information will be available from the quantex CD version 04.
We would like to anticipate these implementations by attaching here the new control values insert sheet of the two precedent lots together with the note.

quantex HbA_{1c}

REF 3000-2314

Note:

The **quantex HbA_{1c}** were correlated to an NGSP-certified High Performance Liquid Chromatography (HPLC) method. The correlation results were:

$$\text{NGSP}_t = 0.776\text{BK} + 0.095 \quad \text{or} \quad \text{BK} = 1.289 \text{NGSP}_t - 0.123$$

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