**biokit**

08186 Llissà d'Amunt  
Barcelona  
Spain

Tel.:+ 34 93 860 90 00  
Fax:+34 93 860 90 17  
e-mail: biokit@biokit.com  
www.biokit.com

Date: 10-10-2007

Ref: 246/MKT/RT10

**SUBJECT MONOGEN Worldwide Presentation**

Up to now, Biokit has offered its own well-known latex test on slide for mononucleosis under two different names:

- **Monogen**, for the US market, available in two different presentations, 50-test kit and 150-test kit (CODE 300021001 and 300021002) and
- **Monolatem** for the rest of the world (ROW) available in a 50-test kit presentation (CODE 3000-1001)

Taking into account that both presentations contain exactly the same reagents, we have decided to use a single name for worldwide use, which will be CE marked and FDA cleared at the same time.

The definitive name will be **Monogen** and the CODE number is maintained as up to now for the ROW version.

The definitive codes and presentations are:

CODE 3000-1001. Monogen. 50-test kit

Containing:

1 x 1.4 ml reagent

1x1 ml positive control

1x1 ml negative control

9 x 6 disposable slides

CODE 3000-1002. Monogen. 150-test kit

Containing:

3 x 1.4 ml reagent

1x2 ml positive control

1x2 ml negative control

27 x 6 disposable slides

Please find attached the new insert and the revised MDSS of this unified reagent:

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---

## monogen

**Rapid latex particle agglutination test on slide for the qualitative and semiquantitative detection of infectious mononucleosis heterophile antibodies in serum or plasma. monogen aids in the diagnosis of infectious mononucleosis.**

---

## Summary

Infectious mononucleosis (IM) is an acute infectious disease of viral etiology. The most frequent symptoms are fever, sore throat, tender lymphadenopathy, anorexia, malaise, headache and myalgia. Splenomegaly occurs in most patients. A macular, maculopapular or petechial rash occurs in up to 50% of the cases, but such rashes occur most commonly in patients who have been treated with ampicillin.

The complications of IM include secondary bacterial pharyngitis, rupture of the spleen, autoimmune hemolytic anemia, autoimmune thrombocytopenia, myocarditis, hepatitis and central nervous system involvement with meningoencephalitis or transverse myelitis. Fatal fulminant IM or acquired hypogammaglobulinemia is rarely seen.

The diagnosis made on clinical history and symptomatology alone is difficult. Numerous cases in which IM has been misidentified with other non-related viral and bacterial diseases have been cited (1). For this reason, hematologic and serologic tests are very helpful in diagnosis.

In 1932, Paul and Bunnell (2) noted that sera from patients with IM have heterophile antibodies to sheep erythrocytes. Also were described agglutinins to red blood cells from other mammals (3,4). The proteins responsible for this agglutination are glycoproteins from red cell membranes called Paul-Bunnell antigen by several authors. Studies made on these glycoproteins show that those purified from bovine red blood cells are the most sensitive to IM heterophile antibodies. Heterophile antibodies to sheep erythrocytes (which are different from those present during IM) may also be detected in sera from normal people, from individuals who have received injections of serum, and others (3,5).

Traditionally the IM heterophile antibodies have been distinguished from other heterophile antibodies by a "differential" absorption test (6,7) with bovine red blood cells and guinea pig kidney tissue. Now, the use of a purified Paul-Bunnell antigen attached to latex particles provides a simple method with improved sensitivity for the specific detection of heterophile antibodies associated with IM.

In 1968, the etiologic agent of IM was described (8). It was called the Epstein-Barr virus (EBV), a member of the herpes virus group. Subsequently, several serologic techniques involving EBV-related antigens have been developed.

The mode of transmission of IM appears to be intimate salivary contact, salivary contamination of eating and drinking vessels and airborne dissemination of EBV (9).

## Principle

The **monogen** reagent is a suspension of polystyrene latex particles of uniform size coated with highly purified Paul-Bunnell antigen from bovine red cell membranes. The degree of purity of the antigen is such that **monogen** only reacts with IM heterophile antibodies. For this reason, "differential" absorptions are not necessary.

Latex particles allow visual observation of the antigen-antibody reaction. If IM heterophile antibodies are present in the sample, the latex suspension changes its uniform appearance and a clear agglutination becomes evident.

## Reagents

- a) **Latex reagent:**  
Suspension of polystyrene latex particles coated with Paul-Bunnell antigen in a buffer. Contains < 0.1% sodium azide.
- b) **Positive control:**  
Rabbit IgG anti-Paul-Bunnell antigen diluted in a buffer. Contains < 0.1% sodium azide.
- c) **Negative control:**  
Non reactive diluted human serum. Contains < 0.1% sodium azide.

## Precautions

**monogen** is intended for IN VITRO diagnostic use.

The reagents contain sodium azide as a preservative. Azides may react with metal plumbing forming explosive components. Upon disposal, please flush abundantly with water.

All human source material used in the preparation of reagents has been tested by an FDA approved method for the presence of HIV 1/2 and HCV antibodies, as well as for HBsAg and found to be negative.

**WARNING: POTENTIALLY BIOHAZARDOUS MATERIAL.**

Because no test method can offer complete assurance of the absence of infectious agents, the reagents should be handled carefully (10).

Dispose all used materials in a suitable biohazardous waste container.

## Storage

The reagents will remain stable through the expiration date shown on the label if stored between 2-8°C. Do not freeze. The reagents can be damaged by improper handling, especially temperature extremes.

The latex reagent, once shaken, must be uniform without visible clumping. When stored a slight sedimentation may occur and should be considered normal. Do not use reagents if they become contaminated.

The reagent dropper dispenses drops of 28 µl ± 10%. The dropper must be held perpendicular to the slide surface and a single drop allowed to fall. Do not use another dropper without previously checking the volume of the drop.

## Available packaging

- **REF** 3000-1001, kit 50 tests.  
1 x 1.4 ml reagent, 1 ml positive control, 1 ml negative control and 9 x 6 disposable slides.
- **REF** 3000-1002, kit 150 tests.  
3 x 1.4 ml reagent, 2 ml positive control, 2 ml negative control and 27 x 6 disposable slides.

## Material required not provided

Normal saline (0.9% NaCl, only for semiquantitative test), automatic pipettes, rotator, stirrers and timer.

## Sample collection

Serum:

Use fresh serum collected by centrifuging clotted blood. If the test cannot be carried out on the same day, it may be stored between 2 and 8°C for no longer than 8 days after collection. For longer periods the samples must be frozen (-20°C).

Plasma:

Collect blood into a tube containing anticoagulant (EDTA). Other anticoagulants should be evaluated before use. Centrifuge to separate plasma from cellular elements. Test the specimen within 24 hours of blood collection.

Do not use hemolyzed or contaminated samples.

## PROCEDURE

Quality control: Before performing a set of determinations it is advisable to test the latex reagent with each of the controls included in the kit, following the steps outlined in the section QUALITATIVE TEST. The reaction between the positive control and the reagent should show a clear agglutination different from the uniform appearance of the negative control. If these results are not obtained, do not use the kit.

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To assure proper delivery the reagent dropper must be held vertically and a single drop allowed to fall.

## QUALITATIVE TEST

- Allow the reagents to reach room temperature (20 - 30°C).
- Gently shake the reagent vial to disperse and suspend the latex particles in the buffer solution. Vigorous shaking should be avoided.
- Place 50 µl of the sample (or one drop of control) onto one section of the slide.
- Shake the reagent vial and add one drop of reagent next to the drop of sample.
- Mix both drops with a stirrer covering the whole surface of the slide section.
- Rotate the slide for 3 minutes manually or on a rotary shaker set at 80-100 rpm.
- Observe for agglutination.

## Interpretation of the results

POSITIVE REACTIONS:

Presence of agglutination. Sample contains IM heterophile antibodies.

- 3+ Large clumping with clear background.
- 2+ Moderate clumping with fluid slightly opaque in background.
- 1+ Small clumping with opaque fluid in background.

NEGATIVE REACTIONS:

Absence of agglutination, uniform suspension.

## SEMIQUANTITATIVE TEST

- Place 50 µl of normal saline onto slide sections 1 through 6.
- Using an automatic pipette, place 50 µl of the sample onto slide section 1.
- Using the same pipette take in and release the sample and the saline solution on section 1 several times until they are well mixed.
- Take 50 µl of the mixture made on section 1 and transfer it to section 2.
- Repeat the aforementioned operations to obtain a thorough mixing of reagents, through section 6, thereafter discarding 50 µl.
- Test each dilution as described in the section QUALITATIVE TEST.

## Interpretation of the results

The IM titer of a sample is the highest dilution showing a positive result.

## Limitations of the procedure

- The results should be interpreted in light of the clinical, hematological and serological information of the patient.
- Occasionally detectable levels of heterophile antibodies are late in developing in patients symptomatic for IM. If symptoms persist it is recommended to repeat the assay in several days. Some patients may remain persistently negative, especially children and adolescent. It has been reported that only 80 to 90% of adults and less than 50% of young children develop heterophile antibodies.
- False positive reactions when testing kits for diagnosis of IM heterophile antibodies have been reported in serum samples collected from patients with recent cytomegalovirus, hepatitis A virus, parvovirus and leptospira infection.
- Detectable levels of heterophile antibodies may persist for months, and more rarely for years, in some individuals.
- Although titers of heterophile antibodies have little relation with the severity of infection, the semiquantitative procedure can be used to follow the evolution of the disease.

## Expected values

Different studies (11,12) of presence of IM heterophile antibodies in blood donors show that the incidence of the disease ranges from 0.9 to 1.7% of population. As presence of heterophile antibodies indicates a relatively recent infection, these results suggest that the true incidence of the disease is higher than the number of diagnosed cases.

## Performance characteristics

Evaluations comparing **monogen** to a commercially available differential hemagglutination test were performed to determine the sensitivity and specificity of the reagent. 106 positive sera and 114 negative sera were selected by a differential hemagglutination slide test for the study. Discrepancies between the results given by **monogen** and the hemagglutination test were resolved using Epstein-Barr Virus (EBV) specific serological assays.

In these assays, the titers of specific antibodies to the EBV capsid antigen (both IgG and IgM), EBV early antigen (both diffused -D- and restricted -R-) and EBV nuclear antigen were determined. The results of these assays specified whether an EBV infection was recent or acute, in which case the serum was considered positive, or if antibodies were absent or relics of an old infection, in which case the serum was considered negative. Twelve of the 100 nondiscrepant positive sera and five of the 106 nondiscrepant negative sera were also analyzed using these same EBV specific serological assays. In all cases the EBV specific assay confirmed the positivity or negativity of the samples.

Compared with hemagglutination, **monogen** was found to have a sensitivity of 94% and a specificity of 93%. Assuming that the concordant results of the EBV serology performed on nondiscrepant sera apply to all the samples tested, it can be inferred that the sensitivity of **monogen** relative to EBV specific tests is 99% and its specificity relative to the same is 93%. There was only one result negative to **monogen** and positive to EBV serology. It was determined to be a recent infection and not an acute infection due to the absence of anti-VCA IgM and nuclear antigen. For this reason this negative result cannot be interpreted as prozone effect.

In a separate study **monogen** was compared to a qualitative red horse cell slide test involving a total of 224 EDTA plasma samples. There was complete agreement in test results which included 51 positive and 173 negative samples. Overall, the results of this study indicate clearly that **monogen** is highly sensitive and specific for the diagnosis of IM.

A panel of 10 positive serum samples was tested on three consecutive days using the semiquantitative technique. The results of the study indicate that **monogen** has 100% precision. The error of repeated estimations was expected to be only one doubling dilution.

## References

1. Davidsohn I, Stern K and Kashiwagi C. The Differential Test for Infectious Mononucleosis. *Amer J Clin Pathol* 21: 1101-1113, 1951.
2. Paul JR and Bunnell WW. The Presence of Heterophile Antibodies in Infectious Mononucleosis. *Amer J Med Sci* 183: 90-104, 1932.
3. Beer P. The Heterophile Antibodies in Infectious Mononucleosis and after the Injection of Serum. *J Clin Invest* 15: 591-599, 1936.
4. Bailey GH and Raffel S. Hemolytic Antibodies for Sheep and Ox Erythrocytes in Infectious Mononucleosis. *J Clin Invest* 14: 228-244, 1935.
5. Forssman J. Die Herstellung Hochwertiger Spezifischer Schaffhämolyse ohne Verwendung von Schafblut. *Biochem Ztschr.* 37: 78-115, 1911.
6. Davidsohn I. Serologic Diagnosis of Infectious Mononucleosis. *JAMA* 108: 289-295, 1937.
7. Stuart CA, Burgess AM, Lawson HA and Wellman HE. Some Cytologic and Serologic Aspects of Infectious Mononucleosis. *Arch Intern Med* 54: 199-214, 1934.
8. Henle G, Henle W and Diehl V. Relation of Burkitt's Tumor-Associated Herpes-Type Virus to Infectious Mononucleosis. *Proc Nat Acad Sci USA* 59: 94-101, 1968.
9. Hoagland RJ. The Transmission of Infectious Mononucleosis. *Amer J Med Sci* 229: 262-272, 1955.
10. Biosafety in Microbiological and Biomedical Laboratories. CDC/NIH manual, 4th Edition, 1999
11. Barret AM. The Serological Diagnosis of Glandular Fever (Infectious Mononucleosis): a New Technique. *J Hyg* 41: 330, 1941.
12. Virtanen, S. Incidence of Infectious Mononucleosis Antibodies in Blood Donors. *Acta Pathol Microbiol Immunol Scand* 56: 53-56, 1962.

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### 9. Physical and chemical properties

**Appearance:**

Latex Reagent: white suspension  
 Positive Control: colourless liquid  
 Negative Control: colourless liquid

**Odour:** Odourless

**Boiling point:** N/A

**Melting point:** N/A

**Flash point:** N/A

**Ignition temperature:** N/A

**Explosion limits:** N/A

**Vapour pressure:** N/A

**Density:** N/A

**Viscosity:** N/A

**Solubility in water:** soluble

### 10. Stability and reactivity

**Conditions and materials to avoid:** N/A

**Hazardous reactions:** N/A

**Hazardous decomposition products:** N/A

### 11. Toxicological information

Data not available.

### 12. Ecological information

Data not available.

### 13. Disposal considerations

Remains of samples, reagents and controls should be collected in a recipient for this purpose and autoclaved 1 hour at 121°C. Observe all Federal, State and local laws.

### 14. Transport information

**Special requirements:** None

### 15. Regulatory information

According to 1999/45/EC Directive and 91/155/EEC Directive and following modifications.

BIOKIT S.A. 08186 Llissà d'Amunt (Barcelona) Spain Phone: (34) 93 860 9000 Fax: (34) 93 860 9029 e-mail: biokit@biokit.com	<b>SAFETY DATA SHEET</b>	Date: 10-10-07	Rev.: 01
	Code: 3000-1001 3000-1002		Pag.: 3 of 3

#### 16. Other information

**Intended use:** Rapid latex particle agglutination test on slide for the qualitative and semiquantitative detection of infectious mononucleosis heterophile antibodies in serum or plasma.

For IN VITRO diagnostic use.

All reagents in this kit contain sodium azide as a preservative. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Use proper disposal procedures.

Each donor unit used in the preparation of the negative control of this kit was tested by an approved method for the presence of the antibodies to HIV and HCV as well as for HBsAg and found to be negative. **WARNING:** Because no test method can offer complete assurance that HIV, HCV, HBsAg or other infectious agents are absent, the controls of this kit should be handled carefully.

The above information is believed to be correct, but does not purpose to be all inclusive and should be used only as a guide. BIOKIT S.A. shall not be held for any damage resulting from handling or use of the product.