



08186 Llíssà d'Amunt
 Barcelona
 Spain

Tel.:+ 34 93 860 90 00
 Fax:+ 34 93 860 90 17
 e-mail: biokit@biokit.com
www.biokit.com



Date: 07/01/04

Ref: 149/MKT/Q/53

SUBJECT BQA New Software version

A new main software version for the BQA instrument is already available (2550164125).

This new version is totally compatible with the last one (2550164124) and user interface software version (2550187106) so if you have this last PC software version is not necessary to change it.

The process to update the Main version is exactly the same than other versions, using FFFTP and Telnet software.

The table below is the current summary of all the PC software versions available for BQA together with its interrelationship with the instrument software versions.

	PC	Main	Sub
	2550187	2550164	2550165
05/07/02	102	113	100
21/04/03	104	115	100/101
06/06/03	105	119	100/101
01/08/03	"	"	"
15/10/03	T 106	122	"
01/11/03	106	124	"
26/11/03	"	125	"

If you need to update your BQA instrument or you have any doubt or question, please don't hesitate to contact with Oscar Marquez (omarquez@biokit.com).