



08186 Llíssà d'Amunt
Barcelona
Spain

Tel.:+ 34 93 860 90 00
Fax:+ 34 93 860 90 17
e-mail: biokit@biokit.com
www.biokit.com

Date: 16/02/05

Ref: 178/MKT/E/58

SUBJECT **Best 2000 – Instruments
Serial Number 1DXB0755 -
1DXB0817 - Deep Well Plate**



There was a batch of the Injection Molded Deep Well plate bases that have been identified as being undersized. This can cause the deep well plates to fit tighter than intended.

Most of the warped bases were caught before they were shipped, however a small quantity may have shipped out. If you have problems fitting the deep well plates onto their bases, please contact Biokit for a replacement.