



08186 Lliçà d'Amunt  
Barcelona  
Spain

Tel.:+ 34 93 860 90 00  
Fax:+ 34 93 860 90 17  
e-mail: biokit@biokit.com  
www.biokit.com

**Date:** 27/04/05

**Ref:** 189/MKT/Q/66

**SUBJECT** BQA New packing list



**We would like to inform you that a new BQA parking list has been created:**

- 1 x 5 litres bottle for wash solution has been replaced for 3 x 2litres so there is one for each wash solution input.
- Plastic tube for high concentrated waste has been replaced for silicon tube. Plastic and silicon tubes have both sides correctly labelled.
- New operator's manual.
- New Screw driver.

**This new version will be included on BQA instruments from serial number 5802-0165.**

**See the attached table for further details.**

**Please, don't hesitate to contact us if you need more information.**

## BQA/IL Packing List

Apr-05

No	品名	Part Name	型式・Type	
1	取り扱い説明書	Operator's manual		1
2	電源ケーブル	Power cable for analyzer	25S1031	1
3	10BASETケーブル	LAN cable	NWNMC5E-STN-CSMB-GY-2	1
4	ハロゲンランプ	Halogen lamp	NKF64258 12V 20W	1
5	ヒューズ	Glass tube fuse 6.3A/250V	218 06.3 6.3A/250V	2
6	タイムラグヒューズ	Glass tube fuse 1.6A/250V	2180 1.6 1.6A/250V	1
7	タイムラグヒューズ	Glass tube fuse 5A/250V	218005 5A/250V	3
8	ミゼットヒューズ	Glass tube fuse 10A/125V	217010 10A/125V	1
9	ヒューズ	Glass tube fuse 3.15A/250V	218 3.15 3.15A/250V	1
10	純水用チューブ組品	Plastic tube for purified water	25-012-4014	1
11	高濃度廃液用チューブ組品	Silicon tube for high conc. waste	25-012-4017	1
12	低濃度廃液用チューブ組品	Plastic tube for low conc. waste	25-012-4016	1
13	洗剤-1用チューブ組品	Plastic tube for wash solution-1	25-012-4011	1
14	洗剤-2用チューブ組品	Plastic tube for wash solution-2	25-012-4012	1
15	洗剤-3用チューブ組品	Plastic tube for wash solution-3	25-012-4013	1
16	チップ交換治具	Syringe tip insertion tool		1
17	PTFEチップ	Syringe tip (TEF010)	TEF010	1
18	PTFEチップ	Syringe tip (TEF050)	TEF050	1
19	PTFEチップ (Oリング付き)	Syringe tip (TEF250) with O ring	TEF250	3
20	PTFEチップ (Oリング付き)	Syringe tip (TEF500) with O ring	TEF500	3
21	+ドライバー	Screw driver (#1-3)	No123-C-100	1
22	六角レンチ	Hexagonal wrench (1.5 mm)		1
23	六角レンチ	Hexagonal wrench (0.9 mm)		1
24	洗剤 フルクリーンNo. 10	Wash solution No.10-2		1
25	2 Lタンク	2L Plastic tank	25-012-4031	3
26	10Lタンク	10L Plastic tank		1
27	20Lタンク	20L Plastic tank		2
28	ASPTレイ組品	Sample tray Assy.	25AP-X-4425	1
29	ASP蓋組品	ASP Lid	25AP-X-4423	1
30	RCUTレイ組品	Reagent bottle tray Assy.	25AP-X-4525	1
31	RCU蓋組品	RCU Lid	25AP-X-4524	1