



08186 Lliçà d'Amunt  
Barcelona  
Spain

Tel.:+ 34 93 860 90 00  
Fax:+ 34 93 860 90 17  
e-mail: [biokit@biokit.com](mailto:biokit@biokit.com)  
[www.biokit.com](http://www.biokit.com)

**Date:** 14/02/06

**Ref:** 212/MKT/Q/72

**SUBJECT PRODUCT MODIFICATION**

## Quantex Ferritin standard multipoint

### Informative Notice:

In order to avoid any possibility of error in the quantex Ferritin calibration, we are taking out the 1000 ng/mL calibrator from the calibrator set *quantex FERRITIN standard multipoint* (COD. 3000-2223).

This standard was included up to know but was not used in the calibration.

This will be applied from the lot. A-1506

