



**biokit**

08186 Llissà d'Amunt  
Barcelona  
Spain

Tel.:+ 34 93 860 90 00  
Fax:+34 93 860 90 17  
e-mail: biokit@biokit.com  
www.biokit.com

Date: 28.06.07

Ref: 239/MKT/Q/50

**SUBJECT FROZEN STABILITY CLAIM FOR  
QUANTEX FERRITIN/MYOGLOBIN/IgE  
CONTROL I/II**



With the aim of improving the performance of our reagents, and being more competitive, we have performed all the necessary experiments using different lots of our **Quantex FERRITIN/MYOGLOBIN/IgE CONTROL I/II** (Cod. 3000-2222) to support a new stability claim.

We are happy to announce that, once rehydrated, this control material can be stored fractionated in aliquots and frozen up to **5 months**. The material can be thawed only once.

The experiments were performed using 0.5-ml fractions and cryotubes. It is very important to use cryotubes with caps to keep the sample in good condition, avoiding evaporation and/or contamination.

This new statement data will appear in the next revision of the corresponding insert sheet.