



**biokit**

08186 Llíssà d'Amunt  
Barcelona  
Spain

Tel.:+ 34 93 860 90 00  
Fax:+34 93 860 90 17  
e-mail: biokit@biokit.com  
www.biokit.com

Date: 18/06/07

Ref: 237/MKT/Q/50

**SUBJECT CALIBRATION STABILITY  
QUANTEX RF PLUS**

## Longer Calibration Stability: 30 days on ILab 600/650



Informative Notice:

With the aim of improving the characteristics of our reagents, we have carried out all the necessary experiments using different lots of our **quantex RF plus** reagents to support a new calibration stability claim in the ILab 600/650 series.

Having performed these tasks and achieving positive results, we are happy to announce that the calibration stability period will now be **30 days** instead of the original 7 days.

This new data statement will appear in the next revision of the corresponding insert sheet.

The sentence:

*“Recalibrate every 7 days, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.”* appearing in the current insert version, will now read:

*“Recalibrate every 30 days, when a new lot of reagent is used, when control recovery falls out of the expected range or when adjustments are made to the instrument.”*